

# Adapters

## Monopolar adapters

### Monopolar adapter, VIO, ICC, ACC, Standard (MO 9/5 mm)



for connecting plug Ø 4mm (single-pole)  
adapter for endoscopy, for STORZ, WOLF, OLYMPUS

**No. 20183-026**


= 1

4 kVp

30183-128



### Monopolar adapter, VIO, ICC, ACC, Standard (MO 9/5 mm)



for connecting plug international (3 Pin)

**No. 20183-028**

= 1

4 kVp

30183-128



### Monopolar adapter, VIO, ICC, ACC, Standard (MO 9/5 mm)



for connecting plug Ø 3 mm

**No. 20183-031**

= 1

4 kVp

30183-128



### Monopolar adapter, VIO, ICC, ACC, Standard (MO 9/5 mm)



for connecting plug Ø 8 mm (Bovie Jack)

**No. 20183-034**


= 1

4 kVp

30183-128



### Monopolar adapter, VIO, ICC, ACC, Standard (MO 9/5 mm)



for standard Martin accessories (MO 8/4 mm)

**No. 20183-047**


= 1

4 kVp

30183-128



### Monopolar adapter, VIO, ICC, ACC, Standard (MO 9/5 mm)



for connecting plug 3 x 2,35 mm Pin (Olsen)

**No. 20183-038**

= 1

4 kVp

30183-128



## Monopolar adapters

### Monopolar adapter, VIO, ICC, ACC, Standard (MO 9/5 mm)



for connection of Lamidey handle (5 poles)

**No. 20183-049**



= 1

0.5 kVp



30183-128



### Monopolar adapter, VIO, ICC, ACC, Standard (MO 9/5 mm)



for APC connection (5 poles)

**No. 20183-063**



= 1


4 kVp



30183-135



### Monopolar adapter, VIO, ICC, ACC, International Ø 8 mm (Bovie Jack)



for APC connection (5 poles)

**No. 20183-064**



= 1


4 kVp



30183-135



### Monopolar adapter, T-Series, Standard (MO 10/4 mm)



for connecting plug international (3 Pin)

**No. 20183-029**



= 1

4 kVp



30183-128



### Monopolar adapter, T-Series, Standard (MO 10/4 mm)



for ERBE connecting plug (single-pole) VIO, ICC, ACC, Standard (MO 9/5 mm), without hand activation

**No. 20183-043**



= 1

4 kVp



30183-128



### Monopolar adapter, EDER 100 L (MO 5/2 mm)



for connecting plug Ø 8 mm (Bovie Jack)

**No. 20183-035**



= 1

4 kVp



30183-128



# Adapters

## Bipolar adapters

### Bipolar adapter, VIO, ICC, ACC, T-Series, Standard (BI 8/4 mm)



for connecting plug international (2- Pin 19 mm)

**No. 20183-017**



= 1

0.8 kVp



30183-128



### Bipolar adapter, VIO, ICC, ACC, T-Series, Standard (BI 8/4 mm)



for Aesculap bipolar connecting plug (old version)

**No. 20183-019**



= 1

0.5 kVp



30183-128



### Bipolar adapter, VIO, ICC, ACC, Standard (BI 8/4 mm)



for Martin, Berchtold, Aesculap connecting plug,  
Standard (BI 9/2.5 mm)

**No. 20183-025**



= 1


0.8 kVp



30183-128



### Bipolar adapter, ERBE units, Standard (BI 8/4 mm)



for Valleylab connecting plug (2-Pin 28 mm)

**No. 20183-065**



= 1


0.5 kVp



30183-128



### Bipolar adapter, VIO, ICC, ACC, T-Series, Standard (BI 8/4 mm)



for connecting plug international (2 x 4 mm socket),  
length 30 cm

**No. 20183-066**



= 1

1 kVp


max 95 °C



30183-153

max 138 °C

### Bipolar adapter, ACC, International (2-Pin 19 mm)



for standard ERBE connecting plug (BI 8/4 mm)

**No. 20183-010**



= 1

0.8 kVp




30183-128



## Bipolar adapters

### Bipolar adapter, ERBE units, International (2-Pin 22 mm)



for standard ERBE connecting plug (BI 8/4 mm)

**No. 20183-053**



= 1


0.8 kVp



30183-128



### Bipolar adapter, Valleylab units (2-Pin 28 mm)



for standard ERBE connecting plug (BI 8/4 mm)

**No. 20183-055**



= 1


0.5 kVp



30183-128



### Bipolar adapter, ERBE units, International (2-Pin 22 mm)



for Valleylab connecting plug (2-Pin 28 mm)

**No. 20183-062**



= 1

0.8 kVp



30183-134



## Patient plate adapters

### Patient plate adapter, VIO, ICC, ACC, T-Series, Standard (NE 6 mm)



for international patient plate cable (NE 2-Pin)

**No. 20183-014**



= 1

4 kVp



30183-128



# Adapters

## Ignition test adapters

### APC ignition test adapter for VIO, ICC, Standard (NE 6 mm)



**No. 20183-070**

= 1

max 95 °C

30183-126

max 134 °C

### APC ignition test adapter for VIO, ICC, International (NE 2-Pin)



**No. 20183-071**

= 1

max 95 °C

30183-126

max 138 °C

## Irrigation adapters

### Irrigation adapter for flexible APC probes with Luer-Lock



for APC probes with APC connecting cable 20132-158

**No. 20132-203**

= 5

max 95 °C

30132-209

max 138 °C

### Irrigation adapter for flexible APC probes with Luer-Lock



for APC probes with APC connecting cables 20132-046

**No. 20132-095**

= 5

max 95 °C

30132-209

max 138 °C

### Irrigation for adapter APC applicators with Luer-Lock



for APC applicators with APC handles

**No. 20132-087**

= 1

max 95 °C

30132-209

max 138 °C

### Irrigation connector for micro forceps with Luer-Lock



**No. 20191-404**

= 1

max 95 °C

30191-799

max 138 °C

## Foot switch adapters

Foot switch adapter from VIO D/S to ICC two-pedal foot switch (4 poles), length 0.7 m



for OR-System e.g. DaVinci

**No. 20140-004**

= 1

30140-212



Foot switch adapter from VIO to ICC one-pedal foot switch (7 poles), length 0.7 m



for OR-System e.g. DaVinci

**No. 20140-007**

= 1

30140-212



## Adapters for special purpose

Adapter Bipolar Resection

Only for VIO with MF socket from version 2.2.1 or 1.8.0



**No. 20183-110**

= 1

30183-168



**ELECTROSURGERY**  
**VESSEL SEALING**  
**ARGON PLASMA COAGULATION**  
**CRYO SURGERY**  
**WATERJET SURGERY**

**ERBE Elektromedizin GmbH**  
**Waldhoernlestrasse 17**  
**72072 Tuebingen**  
**Germany**  
**Phone +49 7071 755-0**  
**Fax +49 7071 755-179**  
**info@erbe-med.com**  
**[www.erbe-med.com](http://www.erbe-med.com)**



**ERBE**

*Perfection for Life*